



Certificate in relation to air transportation of Mobility Aids with Lithium Batteries by passengers with walking disabilities

We hereby certify that the product

erivo R10 (and R11 with 11Ah battery)

including the **lithium-ion battery** INR18650-4P7S-7221 with **11Ah** fulfil the requirements of

IATA chapter 2.3.2.4

and therefore the system may be carried as checked baggage without restriction.

Only the battery has to be removed for transportation (protected from unintentionally switching on, damage and from short circuit by other luggage e.g., by placing the battery in a protective pouch) and **carried** in the **passenger cabin**.

The lithium-ion batteries are tested according UN Manual of Tests and Criteria" Section 38.3.

The rated energy of the battery is with 284,9 Wh **below** the IATA limit of **300 Wh** for batteries which can be removed by the user.

The pilot in command must be informed about the location of the battery.

We strongly recommend contacting the airline sufficiently in advance to request approval for transport.

This certificate does not apply to damaged or defect batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31st December 2024**.

Albstadt, 10th of November 2023

Alber GmbH

A handwritten signature in blue ink, appearing to read "Joost Beltman", with a long horizontal flourish extending to the right.

Joost Beltman
Managing Director