



Certificate in relation to air transportation by passengers with walking disabilities

We hereby confirm that the our product

e-motion M15

including the Lithium ion batteries can be carried by passengers with reservation of the airline company under

IATA chapter 2.3.2.4

Condition precedent to meet the above listed requirement is the transportation as checked baggage in a way that an unintentionally switching on and damage by other luggage is avoided.

The batteries must be removed for transportation and carried (protected from damage and from short circuit) in the **passenger cabin**.

The pilot in command has to be informed about the location of the batteries.

The rated energy of the batteries in each wheel is **140,4 Wh**.

The remote control ECS of the e-motion could be transported also under IATA chapter 2.3.5.9 which allows the transportation via cabin baggage.

We strongly recommend contacting the airline sufficiently in advance to request approval for transport.

This certificate does not apply to damaged or defective batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31st December 2023**.

Albstadt, 02nd of December 2022

Alber GmbH

A handwritten signature in blue ink, appearing to read "Joost Beltman", with a long horizontal stroke extending to the right.

Joost Beltman
Managing Director