

Certificate in relation to air transportation by passengers with walking disabilities

We hereby certify that our product

e-motion M25

including the lithium-ion batteries fulfil the requirements of

IATA chapter 2.3.2.4

and therefore are suitable for transportation by passengers according
The pilot in command must be informed about the location of the mobility aid.

The built-in lithium-ion batteries (**150 Wh**) of the e-motion M25, which are not intended
for removal by the user, fulfil the requirements of the tests listed in the UN Manual of
Tests and Criteria" Part III subsection 38.3.

The e-motion M25 can be set into a dedicated "flight mode". In this mode the system
is completely deactivated. While in this mode the e-motion M25 can only be switched
on by a coded procedure. Unintended switching on by luggage is not possible.
Furthermore, the Bluetooth device of the system is completely switched off.

**We strongly recommend contacting the airline sufficiently in advance to
request approval for transport.**

This certificate does not apply to damaged or defective batteries.

Due to the frequent amendments of the respective regulations, which partly occur on
an annual basis this certificate is valid until **31st December 2023**.

Albstadt, 02nd of December 2022

Alber GmbH



Joost Beltman
Managing Director