

## **Certificate in relation to air transportation by passengers with walking disabilities**

We hereby certify that the our product

### **Twion M24 / T24**

including the Lithium ion batteries fulfil the requirements of

### **IATA chapter 2.3.2.4**

and therefore may be carried as checked baggage without restriction.

The built-in lithium-ion batteries of the e-motion M24 /T24, which are not intended for removal by the user, fulfil the requirements of the tests listed in the UN Manual of Tests and Criteria" Part III subsection 38.3.

Condition precedent to meet the above listed requirement is the transportation of the twion wheels in a way that an unintentionally switching on and damage by other luggage is avoided. We recommend using the dedicated Alber transport axle.

**We strongly recommend contacting the airline sufficiently in advance to request approval for transport.**

This certificate does not apply to damaged or defective batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31<sup>st</sup> December 2023**.

Albstadt, 02<sup>nd</sup> of December 2022

Alber GmbH



Joost Beltman  
Managing Director