



One Invacare Way  
Elyria, OH 44035  
Tel: 866 577 3507  
Fax: 216 325 1111

# e-fix<sup>®</sup> ECO

ORDER  QUOTATION

# 1

- Standard
- Optional field
- Mandatory fields:<sup>1</sup>

The order cannot be processed without details!

**Invoice Address** Customer No. .....

(stamp or completely legible) PO No. ....

Company .....

Street .....

City ..... ST ..... ZIP .....

**Shipping Address** (if different from invoice address)

Company .....

Street .....

City ..... ST ..... ZIP .....

## Type of Order

Dealer installation \*Pending service dept approval\*<sup>2)</sup> \$ 0.00  
(Please complete wheelchair selection)

Installation by Alber<sup>3)</sup> \$ 100.00  
\*All rigid rim chairs must be installed at the Alber factory\*

- e-fix ECO (265 lbs max.) \$5,349.00
- Battery pack, 36V/7,65 Ah \$ 0.00

**Mode of Shipment**  Standard

**Order marked for** .....

- Size of drive wheel
- 14" x 2 1/4 flat free tire with solid insert
  - 22" x 1 3/8 flat free tire with solid insert
  - 24" x 1 3/8 flat free tire with solid insert

## Wheelchair Details

Wheelchair manufacturer (e.g. Ki Mobility) .....

Model (e.g. CATALYST 5) .....

Year of construction .....


Variants (e.g. Amputee) .....

Equipped with vehicle retention system  
(e.g. Tiedown or EZ Lock)


Drum brakes<sup>4</sup>  Yes  No  
 Note: Drum brakes must be removed!



**Seat width** (seat tube measured on the outside) ..... inch



**Seat height** (rear, measured to lower edge of seat sling) ..... inch



Standard control unit mounting  Left (instead of right)  Right

## Important Hardware

- Wheelchair mounting bracket \$ 0.00
- 1593206 Anti tipper with jack up function<sup>5)</sup> (pair) \$ 745.00

## Armrest style

- Tubular or swing away
- Desk or full pad style:  With transfer loop  Without transfer loop

**Comments** (mounting, delivery, etc.)

.....

.....

# e-fix<sup>®</sup> ECO

# 2

Standard

Optional field

**Mandatory fields:**<sup>1</sup>

The order cannot be processed without details!

## Control Unit Accessories

1593361	Swing away arm	\$ 736.00
1565836	Protection bumper	\$ 136.00
1566197	Extension for Speed Selection	\$ 101.00
1464165	Attendant mount (for standard controller)	\$ 344.00

### Side:

- Right  
Left

## Joystick Option

<input checked="" type="checkbox"/>	1565576	Standard Joystick	\$ 0.00
<input type="checkbox"/>	1565576	Standard Joystick replacement	\$ 21.00
<input type="checkbox"/>	1565669	Option special, 40mm/1.57 inch	\$ 64.00
<input type="checkbox"/>	1566196	Option soft, 70mm/2.76 inch	\$ 70.00
<input type="checkbox"/>	1591586	Option T-handle	\$ 237.00



Standard



Special



Soft



T-handle

## Intuitive Attendant Control

1566165	Intuitive attendant control unit (Push/Brake aid)	\$ 1,334.00
1467804	Opposing push handle (Recommended for optimal operation)	\$ 133.00
1594174	Intuitive cable extension <sup>6)</sup> (Must be used)	\$ 358.00

Intuitive control mounting hardware (Select style below)

### Side:



### Style:

Right	1489073	Push handle	\$ 382.00
Left	1490671	Stroller handle	\$ 382.00

## External Charging Socket

1594164	Remote Charging Socket <sup>6)</sup>	\$ 213.00
---------	--------------------------------------	-----------

## Spare Battery Pack/Extra Charger

1594565	Battery, 36V/7,65Ah	\$ 876.00
1594566	Battery Bag	\$ 216.00
1593976	Battery Charger	\$ 1,184.00
1519495	Battery Tray (Tilt Chairs)	\$ 619.00

## Wheel Accessories

1566334	Extension for coupling rings	\$ 500.00
---------	------------------------------	-----------

## Spoke Covers for Drive Wheel

1566666	Transparent, 24"	(pair) \$ 171.00
1566203	Transparent, 22"	(pair) \$ 168.00
1566670	Alber Logo, 24"	(pair) \$ 171.00
1566201	Alber Logo, 22"	(pair) \$ 167.00



Anti-Tippers with Jack-Up Function

Contact  
Name: .....

Email: .....

Date .....

Signature of ordering party .....

As at January 1, 2015. All prices are subject to the statutory VAT, and include packaging FOB. The (general Terms and Conditions of Alber USA LLC shall apply; subject to changes.

1) In order to ensure that your order is processed smoothly and promptly, we require ALL the information in the fields shaded in violet  
2) Mounting kits can only be sent to trained specialist dealers (mounting training).  
3) Please send the wheelchair with a copy of this order as soon as possible and carriage-paid to Altec  
4) If drum brakes are installed, these must be removed as they are not compatible with e-fix drive wheels. Standard lever brakes are used instead.  
5) Alber anti-tip support are not compatible with all chairs. Please use chair manufacture anti-tips in these cases.  
6) Must choose either intuitive cable extension OR remote charging socket. Cannot choose both.