



## PRODUCT SELECTION GUIDE

We advise that you take a short moment to assess your needs before ordering your product in order to make the best-possible gain from the mobility benefits. This overview of the most important distinguishing features between the e-motion and twion will help you choose the right product for you.

	e-motion®	twion®
		
<b>Battery capacity / range</b>	approx. 25 km	approx. 15 km
<b>Weight per wheel (incl. battery)</b>	7,8 kg <sup>(1)</sup>	6,3 kg <sup>(1)</sup>
<b>Motor power</b>	Very high (40 Nm torque per wheel)	Medium (20 Nm torque per wheel)
<b>Speed</b>	6 km/h (optional 8,5 km/h) <sup>(2)</sup>	6 km/h (optional 10 km/h) <sup>(2)</sup>
<b>Wheelchair compatibility</b>	Compatible with almost all manual wheelchairs currently on the market	Compatible with almost all manual wheelchairs currently on the market
<b>Max. weight of occupant</b>	150 kg	120 kg
<b>Characteristics and fields of application</b>	<ul style="list-style-type: none"> <li>- Little or unequal arm strength (disbalance)</li> <li>- Limited hand and finger function</li> <li>- Limited upper body stability</li>   <li>- Easily cope with longer distances and inclines</li> <li>- Mobility throughout the day without exhaustion</li> <li>- Therapeutic benefit (training without overload)</li> </ul>	<ul style="list-style-type: none"> <li>- Good coordination and strength of arms</li> <li>- No limitations of hand and finger functions</li> <li>- Good upper body stability</li>   <li>- Expand radius of action</li> <li>- Drive faster</li> <li>- Protect muscles, joints and tendons from permanent strain and wear</li> </ul>
<b>Special usage aids</b>	<ul style="list-style-type: none"> <li>- Adjustable sensitivity of the sensors on the push rims (left and right)</li> <li>- Four pre-defined driving profiles</li> <li>- Two assistance levels per driving profile</li> <li>- Rollback delay allows the safe handling of slopes and ramps</li> <li>- Individual programming of the driving parameters</li> <li>- Cruise Mode function <sup>(2)</sup></li> </ul>	<ul style="list-style-type: none"> <li>- Sport and Eco mode <sup>(2)</sup></li> <li>- Cruise mode function <sup>(2)</sup></li> </ul>

1) Depending on options and accessories (tyres, push-rim etc.).

2) Optionally available with Mobility Plus Package.